

## Anti-Vesicle-associated protein 1-1 antibody

Catalog: PHY2211A

## **Product Information**

**Description:** Rabbit polyclonal antibody

**Background:** VAP27 has high homology to the VAP33 family of SNARE-like proteins from

animals. It may be involved in vesicular transport to or from the ER.

Synonyms: VAP, VAMP/SYNAPTOBREVIN-ASSOCIATED PROTEIN 27-1, ATVAP, VAP27,

VAP27-1, VESICLE ASSOCIATED PROTEIN

Immunogen: KLH-conjugated synthetic peptide (13 aa from C terminal section) derived from

Arabidopsis thaliana VAP (AT3G60600).

Form: Lyophilized

**Quantity**: 150 μg

Purification: Immunogen Affinity Purified

**Reconstitution:** Reconstitution with 150  $\mu$ I of sterile 1 $\times$ PBS (PH=7.4).

"Note: please spin tube briefly prior to opening it to avoid any losses that might

occur from lyophilized material adhering to the cap or sides of the tube".

**Stability &**Use a manual defrost freezer and avoid repeated freeze-thaw cycles.

Storage: 12 months from date of receipt, -20 to -70 °C as supplied.

6 months, -20 to -70 °C under sterile conditions after reconstitution.

1 month, 2 to 8℃ under sterile conditions after reconstitution.

**Shipping:** The product is shipped at 4°C. Upon receipt, store it immediately at the

temperature recommended above.

## **Application Information**

Recommended Dilution: Western Blot (1:1000-1:2000)

Note: Optimal dilutions/concentrations should be determined by the

end user.


**Expected / apparent MW:** 28 kDa

**Predicted Reactivity:** For more species homologues information, please contact tech

support at tech@phytoab.com.



## **Application Example**



Recom: 2.5 ng, 10 ng and 25 ng recombinant protein containing the peptide for immunization and having a molecular mass of 50 kDa.

Electrophoresis: 12% SDS-PAGE

**Transfer:** blotting to NC (nitrocellulose) membrane for 1 h.

**Blocking:** 5% skim milk at RT or 4°C for 1 h.

**Primary antibody:** 1:1000 dilution overnight at 4℃.

Secondary antibody: 1:10000 dilution using Goat Anti-Rabbit IgG H&L (HRP)

(Cat# PHY6000).

Detection: using chemiluminescence substrate and image were captured with

CCD camera.