

Anti-Peroxisome biogenesis protein 12 antibody

Catalog: PHY7126S

Product Information

Description: Rabbit polyclonal antibody

Background: PEX12 is RING finger protein involved in peroxisome biogenesis. It is involved

in peroxisomal import of nitric oxide synthase.

Synonyms: PEX12, ABERRANT PEROXISOME MORPHOLOGY 4, APM4, ATPEX12,

PEROXIN-12

Immunogen: KLH-conjugated synthetic peptide (16 aa from N terminal section) derived from

Arabidopsis thaliana PEX12 (AT3G04460).

Form: Lyophilized

Quantity:150 μgPurification:Serum

Peptide affinity form antibody available upon request at info@phytoab.com.

Reconstitution: Reconstitution with 150 µl of sterile water.

"Note: please spin tube briefly prior to opening it to avoid any losses that might

occur from lyophilized material adhering to the cap or sides of the tube".

Stability &Use a manual defrost freezer and avoid repeated freeze-thaw cycles.

Storage: 12 months from date of receipt, -20 to -70° C as supplied.

6 months, -20 to -70°C under sterile conditions after reconstitution.

1 month, 2 to 8°C under sterile conditions after reconstitution.

Shipping: The product is shipped at 4°C. Upon receipt, store it immediately at the

temperature recommended above.

Application Information

Recommended Dilution: Western Blot (1:1000-1:2000)

Note: Optimal dilutions/concentrations should be determined by the

end user.

Expected / apparent MW: 44 kDa

Predicted Reactivity: Among species analyzed, the sequence of the synthetic peptide used

for immunization is 100% homologues with the sequence in *Solanum tuberosum*, *Nicotiana tabacum*, *Solanum lycopersicum*, and 80-99%


homologues with the sequence in Triticum aestivum, Hordeum



vulgare, Sorghum bicolor, Panicum virgatum, Oryza sativa, Zea mays, Setaria viridis, Brassica napus, Brassica rapa, Populus trichocarpa, Cucumis sativus, Vitis vinifera, Medicago truncatula, Physcomitrium patens, Spinacia oleracea, Glycine max, Gossypium raimondii.

For more species homologues information, please contact tech support at tech@phytoab.com.

Application Example



Recom: 2.5 ng, 10 ng and 25 ng recombinant protein containing the peptide

for immunization and having a molecular mass of 41 kDa.

Electrophoresis: 12% SDS-PAGE

Transfer: blotting to NC (nitrocellulose) membrane for 1 h.

Blocking: 5% skim milk at RT or 4° C for 1 h.

Primary antibody: 1:1000 dilution overnight at 4℃.

Secondary antibody: 1:10000 dilution using Goat Anti-Rabbit IgG H&L (HRP)

(Cat# PHY6000).

Detection: using chemiluminescence substrate and image were captured with

CCD camera.