

Anti-Wheat dwarf virus (WDV) Capsid protein antibody

Catalog: NAV1002-002S

Quantity: 200 µL

Immunogen Information:

Background

Wheat dwarf virus (WDV)

Immunogen

KLH-conjugated synthetic peptide (15 aa from Central section) derived from Wheat dwarf virus WDV Capsid protein (Uniprot: A6GVG4 NCBI: ACP28900). We also have antibodies for different epitopes from the Capsid protein. Please request at info@nanodiaincs.com or https://www.nanodiaincs.com.

Basic Information:

Purification: Serum

Peptide affinity form antibody available upon request at info@nanodiaincs.com.

Clonality: Polyclonal Expected MW: 29 kDa Host: Rabbit

Product Information:

Form: Lyophilized Reconstitution

Reconstitution with 200 µL of sterile water.

"Note: please spin tube briefly prior to opening it to avoid any losses that might occur from lyophilized material adhering to the cap or sides of the tube".

Stability & Storage

Use a manual defrost freezer and avoid repeated freeze-thaw cycles.

12 months from date of receipt, -20 to -70 °C as supplied.

6 months, -20 to -70 °C under sterile conditions after reconstitution.

1 month, 2 to 8°C under sterile conditions after reconstitution.

Shipping

The product is shipped at 4°C. Upon receipt, store it immediately at the temperature recommended above.

Applications Information:

Recommended Dilution: WB (1:1000-1:2000)

Predicted Reactivity: Among species analyzed, the sequence of the synthetic

peptide used for immunization is 100% homologues with


the sequence in Barley dwarf virus.

For more species homologues information, please contact

tech support at info@nanodiaincs.com.



Example



Recom: 2.5 ng, 10 ng and 25 ng recombinant protein containing the peptide for immunization and having a molecular mass of 43 kDa.

Electrophoresis: 12% SDS-PAGE

Transfer: blotting to NC (nitrocellulose) membrane for 1 h.

Blocking: 5% skim milk at RT or 4° C for 1 h.

Primary antibody: 1:1000 dilution overnight at 4℃.

Secondary antibody: 1:10000 dilution using Goat Anti-Rabbit

IgG H&L (HRP) (Cat# PHY6000).

Detection: using chemiluminescence substrate and image

were captured with CCD camera.