

Anti-NADH dehydrogenase [ubiquinone] iron-sulfur protein 3 antibody

Catalog: PHY0516S

Product Information

Description:	Rabbit polyclonal antibody
Background:	Complex I is the largest protein complex of the oxidative phosphorylation system in mitochondrial and it catalyzes NADH-quinone oxidoreduction. Complex I represents the main entrance site for electrons into the respiratory electron transfer chain. In Arabidopsis, Complex I have at least 49 subunits and NAD9 (ATMG00070) is one of the subunit.
Synonyms:	NAD9, NADH DEHYDROGENASE SUBUNIT 9
Immunogen:	KLH-conjugated synthetic peptide (16 aa from C terminal section) derived from <i>Arabidopsis thaliana</i> NAD9 (ATMG00070).
Form:	Lyophilized
Quantity:	150 µg
Purification:	Serum Peptide affinity form antibody available upon request at info@phytoab.com .
Reconstitution:	Reconstitution with 150 µl of sterile water. "Note: please spin tube briefly prior to opening it to avoid any losses that might occur from lyophilized material adhering to the cap or sides of the tube".
Stability &	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.
Storage:	12 months from date of receipt, -20 to -70°C as supplied. 6 months, -20 to -70°C under sterile conditions after reconstitution. 1 month, 2 to 8°C under sterile conditions after reconstitution.
Shipping:	The product is shipped at 4°C. Upon receipt, store it immediately at the temperature recommended above.

Application Information

Recommended Dilution:	Western Blot (1:1000-1:2000) Note: Optimal dilutions/concentrations should be determined by the end user.
Expected/apparent MW:	23 kDa

Research Use Only

Predicted Reactivity:

Among species analyzed, the sequence of the synthetic peptide used for immunization is 100% homologous with the sequence in *Gossypium raimondii*, *Solanum tuberosum*, *Brassica napus*, *Triticum aestivum*, *Glycine max*, *Oryza sativa*, *Vitis vinifera*, *Medicago truncatula*, *Hordeum vulgare*, *Setaria viridis*, *Zea mays*, *Solanum lycopersicum*, *Cucumis sativus*, *Brassica rapa*, *Nicotiana tabacum*, and 80-99% homologous with the sequence in *Spinacia oleracea*, *Populus trichocarpa*, *Chlamydomonas reinhardtii*.

For more species homologues information, please contact tech support at tech@phytoab.com.