

# Anti-Cherry rasp leaf virus (CRLV) Coat protein-VP24, C-terminal antibody

Catalog: NAV1008-001S

Quantity: 200 µL

## **Immunogen Information:**

#### **Background**

Cherry rasp leaf virus (CRLV)

## **Immunogen**

KLH-conjugated synthetic peptide (15 aa from C terminal section) derived from Cherry rasp leaf virus CRLV Coat protein-VP24 (NCBI: YP\_081448). We also have antibodies for different epitopes from the Capsid protein. Please request at info@nanodiaincs.com or https://www.nanodiaincs.com.

#### **Basic Information:**

Purification: Serum

Peptide affinity form antibody available upon request at <a href="mailto:info@nanodiaincs.com">info@nanodiaincs.com</a>.

**Clonality:** Polyclonal Expected MW: 22 kDa Host: Rabbit

#### **Product Information:**

Form: Lyophilized Reconstitution

Reconstitution with 200 µL of sterile water.

"Note: please spin tube briefly prior to opening it to avoid any losses that might occur from lyophilized material adhering to the cap or sides of the tube".

#### Stability & Storage

Use a manual defrost freezer and avoid repeated freeze-thaw cycles.

12 months from date of receipt, -20 to -70 °C as supplied.

6 months, -20 to -70 °C under sterile conditions after reconstitution.

1 month, 2 to 8℃ under sterile conditions after reconstitution.

## Shipping

The product is shipped at 4°C. Upon receipt, store it immediately at the temperature recommended above.

## **Applications Information:**


Recommended Dilution: WB (1:1000-1:2000)

**Predicted Reactivity:** For more species homologues information, please contact

tech support at info@nanodiaincs.com.



## **Example**



Recom: 2.5 ng, 10 ng and 25 ng recombinant protein containing the peptide for immunization and having a molecular mass of 52 kDa.

Electrophoresis: 12% SDS-PAGE

Transfer: blotting to NC (nitrocellulose) membrane for 1 h.

**Blocking:** 5% skim milk at RT or 4℃ for 1 h.

**Primary antibody:** 1:1000 dilution overnight at 4°C.

Secondary antibody: 1:10000 dilution using Goat Anti-Rabbit

IgG H&L (HRP) (Cat# PHY6000).

Detection: using chemiluminescence substrate and image

were captured with CCD camera.